


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Dodecanol - reagent grade, 98%

Product Number: 126799
CAS Number: 112-53-8
MDL: MFCD00004753
Formula: C₁₂H₂₆O
Formula Weight: 186.33 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White (Solid or Crystalline Solid) or Colorless (Liquid)	
Appearance (Form)	Conforms to Requirements
Solid or Crystalline Solid or Liquid	
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 97.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
c = 100mg/ml, CHCl ₃	
Solubility (Color)	Colorless
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000019240

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.