

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Reactive Red 120–Agarose – saline suspension, Type 3000–CL

Product Number: **R0503**

MDL: MFCD00165711

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	Red to Very Dark Red
Appearance (Form)	Suspension
Matrix	Conforms
Cross-linked 4% beaded agarose	
mg/ml binding capacity (BSA)	≥ 5
Ligand Immobilized	≥ 3
micromoles Reactive Red 120 per ml packed gel	
ml suspension/ml packed gel	Conforms
2 ml suspension yields 1 mL packed gel	

Specification: PRD.3.ZQ5.10000007322

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.