

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl 4-hydroxybenzoate - BioReagent, suitable for insect cell culture

Product Number:

H3647

CAS Number:

99-76-3

MDL:


MFCD00002352

Formula:

C₈H₈O₃

Formula Weight:

152.15 g/mol



TEST

Specification

Appearance (Color)

Colorless or White

Appearance (Form)

Powder or Crystals

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

50 mg/mL, EtOH

1H NMR Spectrum

Conforms to Structure

Purity (HPLC)

≥ 99.0 %

Recommended Retest Period

5 Years

Insect Cell Test

Pass

Specification: PRD.2.ZQ5.10000016347

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.