

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Anti-Lectin, Ulex europaeus-I antibody produced in rabbit - whole antiserum, delipidized

Product Number:

U4754

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Form)

Liquid

Appearance (Turbidity)

Clear to Slightly Hazy

Appearance (Color)

Very Light Brown to Very Dark Brown and  
Dark Brown-Red to Very Dark Brown-Red

Protein by Biuret

40 - 100 mg/ml

DID Antibody Titer

1:8-1:16

IHC - Tissue

Human Tonsil

IHC - Fixation

FFPE-Formalin Fixed Paraffin Embedded

IHC - Titer

1:500-1:1000

Recommended Retest Period

-----

5 Years

Specification: PRE.0.ZQ5.10000017851

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.