

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium cyanoborohydride solution – 5.0 M in 1 M NaOH

Product Number:

296945



MDL:

MFCD00003516

Formula:

CH₃BNNa

Formula Weight:

62.84 g/mol

TEST

Specification

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Titration with Na₂S₂O₃

4.75 - 5.50 M

Specification: PRD.0.ZQ5.10000052704

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.