

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum on carbon - Evonik Noblyst® P2057 3% Pt

Product Number:

908169

Pt

TEST	Specification
------	---------------

Appearance (Color)	Black
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Platinum Component	
Loss on Drying	45.0 - 65.0 %
Calculated to dry basis	
Platinum (Pt)	2.7 - 3.3 % wt
Calculated to dry basis	
Activity	
Crotonic acid hydrogenation ml/min.g	
Specification: Report Only	
Product of Supplier	Confirmed
Evonik	
Registered Trademark	Note
Noblyst is a registered trademark of Evonik Degussa GmbH	
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000100621

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.