

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Nutrient broth – Microbiologically tested, NutriSelect® Basic

Product Number: **N7519**

TEST	Specification
Appearance Beige	Conforms
Solubility Clear, Pale to Light Yellow, 0.8% aqueous solution	Conforms
pH c= 0.8% in H ₂ O	6.6 - 7.0
Microbiological Test	Pass
Expiration Dated Product	-----
Loss on Drying	≤ 5 %

Specification: PRD.3.ZQ5.10000035339

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.