

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Lectin, Triticum vulgaris antibody produced in rabbit - fractionated antiserum, lyophilized powder

Product Number:

T4144

TEST**Specification**

Appearance (Form)	Lyophilized Powder
Appearance (Color)	White to Off White
Solubility (Solvent)	Water
Solubility (Volume)	2.0 - 2.0 ml
Solubility (Color)	Colorless to Faint Yellow and Colorless to Faint Brown
Solubility (Turbidity)	Clear to Slightly Hazy
Protein by Biuret	10 - 100 mg/ml
IEP vs Triticum Vulgaris Lect.	Positive
ELISA - Titer	1:50,000-1:75,000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000020936

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.