

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Potassium 2-chloro-4-(trifluoromethyl)pyridine-3-trifluoroborate - 96%

Product Number: 790400

Formula: C<sub>6</sub>H<sub>2</sub>BClF<sub>6</sub>KN



TEST	Specification
------	---------------

Appearance (Color)	White to Faint Brown and Faint Beige and Faint Yellow
Appearance (Form)	Solid or Powder or Crystals or Chunks
Infrared Spectrum	Conforms to Structure
Boron NMR Spectrum	Conforms to Structure
Fluorine NMR Spectrum	Conforms to Structure
Carbon	23.9 - 26.2 %
Purity (HPLC)	≥ 95.5 %

Specification: PRD.0.ZQ5.10000045352

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.