


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Piperonal - ≥99%, FCC, FG

Product Number: W291102
CAS Number: 120-57-0
MDL: MFCD00005828
Formula: C₈H₆O₃
Formula Weight: 150.13 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Purity (GC)	≥ 99.0 %
Solidification Point	≥ 35 °C
Solubility (Turbidity) 1ml/4ml, 70% Ethanol	Clear
Solubility (Color)	Colorless to Faint Yellow
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Organoleptics	Pass
Expiration Date Period 5 Years	-----

Specification: PRD.4.ZQ5.10000003680

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

