


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Methyl trans-cinnamate - ≥98%, FCC, FG

Product Number: **W269808**
CAS Number: 1754-62-7
MDL: MFCD00008458
Formula: C₁₀H₁₀O₂
Formula Weight: 162.19 g/mol



TEST	Specification
Appearance (Color)	White to Pale Yellow
Appearance (Form)	Solid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.00 %
Acid Value	≤ 2.0 ml
Chlorine Compounds	Pass
Solubility (Turbidity) 1g/4ml 80% Ethanol	Clear to Slightly Hazy
Solubility (Color)	Colorless to Faint Yellow
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period 5 Years	-----

Specification: PRD.1.ZQ5.10000021257

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.