


Product Specification

Product Name : Poly(acrylic acid)
Product Number: 323667
CAS Number: 9003-01-4
MDL Number : MFCD00084394
Storage Temperature : +2°C to +8°C



Test	Specification
Appearance (Form)	Powder
Appearance (Color)	White
IR	Conforms to Structure
Brookfield Viscosity 25 deg C	≤ 2000cps

Sigma-Aldrich warrants that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Version Number: PRP.00.ZQ5.1162775

Page 1 of 1

The branding on the header and/or footer of this document may temporarily not visually match the product purchased as we transition our branding. However, all of the information in the document regarding the product remains unchanged and matches the product ordered. For further information please contact mlsbranding@sial.com

