

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium methoxide solution - ACS reagent, 0.5 M CH₃ONa in methanol (0.5N)

Product Number:

403067

CH₃ONa

MDL:

MFCD00012179

Formula:

CH₃NaO

Formula Weight:

54.02 g/mol

| TEST | Specification |
|------|---------------|
|------|---------------|

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Titration with HCl

0.48 - 0.52 M

Turbidity Measurement

Pass

Clarity of Solution Test

Specification: PRD.0.ZQ5.10000048705

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.