

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

**Product Specification**

Product Name:  
Anti-Amyloid Precursor Protein, C-Terminal antibody produced in rabbit – IgG fraction of antiserum, buffered aqueous solution

Product Number:                      **A8717**

Storage Temperature:                -20 °C

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	9.0 - 15.0 mg/ml
WB-Cell Line/Tissue Extract	Rat Brain S1 Frac.
WB Titer	1:4,000-1:8,000
IHC - Tissue	Human Alzheimer's Disease Brain
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:1000-1:2000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000016381

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.