

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification


Product Name:

Triphenylphosphine - ReagentPlus*, 99%

Product Number: T84409 CAS Number: 603-35-0

MDL: MFCD00003043

Formula: C18H15P
Formula Weight: 262.29 g/mol



TEST Specification

Appearance (Color)

Appearance (Form)

Crystals, Flakes or Pellets

Infrared spectrum

Purity (GC)

Carbon
Solubility (Turbidity)

c = 50 mg/ml, CHCl3

Solubility (Color)

White

Conforms to Requirements

Conforms to Structure

≥ 98.5 %

81.2 - 83.7 %

Clear to Very Slightly Hazy

Colorless

Specification: PRD.0.ZQ5.10000023579

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.