


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isoamyltriphenylphosphonium bromide - 98%

Product Number: 347779
CAS Number: 28322-40-9
MDL: MFCD00031594
Formula: C₂₃H₂₈BrP
Formula Weight: 415.35 g/mol



TEST

Specification

Appearance (Color)

Appearance (Form)

Infrared spectrum

Water (by Karl Fischer)

Titration by AgNO₃

White to Light Orange

Crystals or Powder or Granules or

Chunks or Crystals with Chunks

Conforms to Structure

≤ 0.1 %

≥ 97.5 %

Specification: PRD.0.ZQ5.10000034985

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.