

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium chloride – plant cell culture tested, ≥99%

Product Number:

P8041**KCl**

CAS Number:

7447-40-7

MDL:

MFCD00011360

Formula:

ClK

Formula Weight:

74.55 g/mol

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Powder

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

20 g + 150 mL, H₂OTitration with AgNO₃

≥ 99 %

Plant Cell Culture Test

Pass

Specification: PRD.1.ZQ5.10000013982

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.