


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Methanol solution - (Methanol:Dimethyl sulfoxide 1:1 (v/v))

Product Number: **650188**



MDL: MFCD00004595

Formula: CH<sub>4</sub>O

Formula Weight: 32.04 g/mol

TEST	Specification
------	---------------

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Mass Balance	49.0 - 51.0 %
% V/V DIMETHYL SULFOXIDE	
Expiration Date Period	-----
2 years	

Specification: PRD.0.ZQ5.10000027216

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.