

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium sulfide nonahydrate - ≥98.0%

Product Number:**S2006**

CAS Number:

1313-84-4

MDL:

MFCD00149184

Formula:

Na₂S · 9H₂O

Formula Weight:

240.18 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

White to Off White

Appearance (Form)

Powder or Chunks

Solubility (Turbidity)

Clear

at 5 g plus 25 ml of H₂O

Solubility (Color)

Colorless to Very Faint Yellow

Titration

≥ 98.0 %

Purity by redox titration

Specification: PRD.0.ZQ5.10000020179

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.