

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1,4-Diazabicyclo[2.2.2]octane - ReagentPlus[®], ≥99%

Product Number:

D27802

CAS Number:

280-57-9

Formula:

C₆H₁₂N₂

Formula Weight:

112.17 g/mol



TEST

Specification

Appearance (Colour)

White

Appearance (Form)

Powder or Crystals or Granules or Chunks

Infrared Spectrum

Conforms to Structure

Purity (GC)

≥ 99.0 %

Not Corrected for Water Content

Water (by Karl Fischer)

≤ 1.0 %

Specification: PRE.1.ZQ5.10000014398

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.