

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Bovine Collagen Solution - Type I, 6 mg/mL, ≥95%, sterile filtered, BSE-Free, suitable for biomedical research

Product Number: **804622**

Storage Temperature: 2 - 8 °C

| TEST                                 | Specification   |
|--------------------------------------|-----------------|
| Concentration (mg/ml)                | 6.0 - 8.2 mg/ml |
| pH                                   | 1.8 - 2.2       |
| Electrophoretic Pattern              | Conforms        |
| > = 95% within alpha, beta and gamma |                 |
| < = 5% faster than alpha             |                 |
| Endotoxin Level                      | Conforms        |
| < 1.0 EU/ml                          |                 |
| Expiration Dated Product             | -----           |

Specification: PRD.2.ZQ5.10000084260

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.