


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5,6-Dihydro-4-hydroxy-6-methyl-2H-pyran-2-one

Product Number: 309680
CAS Number: 33177-29-6
MDL: MFCD00012356
Formula: 128.13
Formula Weight: 128.13 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White to Beige
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Proton NMR spectra	Conforms to Structure
Purity (HPLC)	≥ 88.0 %
Appearance (Color)	White to Beige
Appearance (Form)	Powder
Proton NMR spectra	Conforms to Structure

Specification: PRD.0.ZQ5.10000018439

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.