


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Chlorobenzyl chloride - 95%

Product Number: 111961
CAS Number: 104-83-6
MDL: MFCD00000914
Formula: C₇H₆Cl₂
Formula Weight: 161.03 g/mol



TEST	Specification
Appearance (Color) White (Solid) or Colorless (Liquid)	White or Colorless
Appearance (Form) Solid or Liquid	Solid or Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 94.0 %

Specification: PRD.0.ZQ5.10000019722

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.