

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Nafion® perfluorinated membrane – Nafion® 115, thickness 0.005 in.

Product Number: **541346**
CAS Number: 31175-20-9
MDL: MFCD00132790



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form) Sheets	Conforms
Miscellaneous Supplier Data Visual Quality	Pass
Miscellaneous Supplier Data Complete Hydrolysis	Pass
Size Thickness: Mils	
Miscellaneous Assay Weight Caliper: g/sq. in	0.130 - 0.190
Miscellaneous Assay Expansion MD	
Miscellaneous Assay Expansion TD	
Vendor Information Nafion is a registered trademark of The Chemours Company FC, LLC	Conforms

Specification: PRD.0.ZQ5.10000067260

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.