


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
6-Aminoquinoline - 98%

Product Number: 275581
CAS Number: 580-15-4
MDL: MFCD00006803
Formula: C₉H₈N₂
Formula Weight: 144.17 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Yellow to Orange to Brown to Green to Grey-Green
Appearance (Form)

Powder or Crystals or Crystalline
Chunk(s) or Granule or Flakes or Solid

Infrared spectrum

Conforms to Structure

Purity (HPLC)

≥ 97.5 %

Purity (Titration by HCL04)

97.5 - 102.5 %

Specification: PRD.1.ZQ5.10000014532

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.