


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
7,8-Dihydroxy-4-phenylcoumarin

Product Number: 576441
CAS Number: 842-01-3
MDL: MFCD02058070
Formula: C₁₅H₁₀O₄
Formula Weight: 254.24 g/mol



TEST

Specification

Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	
GC (area %)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000055426

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.