


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(2-Naphthylmethyl)zinc bromide solution - 0.5 M in THF

Product Number: 499021
CAS Number: 152329-44-7
MDL: MFCD01311438
Formula: C₁₁H₉BrZn
Formula Weight: 286.48 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Purity Titration by EDTA

0.45 - 0.65 M

Specification: PRD.2.ZQ5.10000024980

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.