


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Phenyl-3-pyrazolidinone - 97%

Product Number: 127914
CAS Number: 92-43-3
MDL: MFCD00003094
Formula: C₉H₁₀N₂O
Formula Weight: 162.19 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Off-White to Yellow	
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	64.3 - 69.0 %
Nitrogen	16.8 - 17.8 %
Purity (HPLC)	≥ 96.5 %

Specification: PRD.1.ZQ5.10000005718

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.