


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Scandium(III) oxide - 99.99% trace rare earth metals basis

Product Number: 751863
CAS Number: 12060-08-1
Formula: O3Sc2
Formula Weight: 137.91 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Complexometric EDTA % Sc	63.9 - 66.5 %
ICP Major Analysis Confirms Scandium Component	Confirmed
Purity 99.99% Based On Rare Earth Analysis	Meets Requirements
Trace Rare Earth Analysis	≤ 150.0 ppm

Specification: PRD.1.ZQ5.10000043617

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.