

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Silver – Nanowires, diam.  $\times$  L 25 nm ( $\pm$ 3 nm)  $\times$  15  $\mu$ m ( $\pm$ 2  $\mu$ m), 5 mg/mL (in IPA)

Product Number:

806846

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	Grey
Appearance (Form)	Suspension
ICP Major Analysis	Confirmed
Confirms Silver Component	
Concentration	Confirmed
5 mg/mL	
Diameter	Confirmed
25+ /-3 nm	
Length	Confirmed
15+ /-2 $\mu$ m	
Isopropyl Alcohol	Confirmed
Recommended Retest Period	-----
1 Year	

Specification: PRD.0.ZQ5.10000102803

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

