


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3,4-Dimethoxybenzaldehyde - 99%

Product Number: 143758
CAS Number: 120-14-9
MDL: MFCD00003363
Formula: C₉H₁₀O₃
Formula Weight: 166.17 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Light Yellow or Tan

Appearance (Form)

Powder or Flake or Crystals or
Crystalline Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 98.5 %

Specification: PRD.1.ZQ5.10000013376

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.