

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-PSD95 antibody produced in mouse - clone 7E3-1B8, purified immunoglobulin, buffered aqueous solution

Product Number: **P246**

MDL: MFCD01324382

Storage Temperature: -20 °C

TEST	Specification
------	---------------

Formulation Solution in PBS containing 0.05% Sodium Azide	Conforms
Concentration 1 mg/mL	Conforms
Working Dilution by Western Blot 1:2,000	Conforms
Working Dilution by Immuno-Fluorescence 1:200	Conforms

Specification: PRD.0.ZQ5.10000059938

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.