

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phenethyl propionate ≥98%, FG

Product Number:

W286702

CAS Number:

122-70-3

MDL:


MFCD00027008

Formula:

C₁₁H₁₄O₂

Formula Weight:

178.23 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 °C	1.491- 1.495
Infrared spectrum	Conforms to Structure
Purity (GC)	≥98.00 %
Color test (APHA)	APHA < 100
Arsenic (As)	≤3 ppm
Cadmium (Cd)	≤1 ppm
Mercury (Hg)	≤1 ppm
Lead (Pb)	≤10 ppm
Expiration Date Period 5 Years	-----

Remarks:

Specification Date : 06/28/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.