


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1,2-Bis[(2S,5S)-2,5-dimethylphospholano]ethane(cyclooctadiene)rhodium(I) tetrafluoroborate

Product Number: **675547**
CAS Number: 213343-65-8
MDL: MFCD04038742
Formula: C₂₂H₄₀BF₄P₂Rh
Formula Weight: 556.21 g/mol



TEST	Specification
Appearance (Color)	Orange
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Proton NMR spectrum	Conforms to Structure
Phosphorus NMR Spectra	Conforms to Structure
Carbon	45.1 - 49.9 %
ICP: Confirms P, Rh and B Comp	Confirmed

Specification: PRD.0.ZQ5.10000054129

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.