

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Lithium hexafluoroarsenate(V) - 98%

Product Number: 308315
CAS Number: 29935-35-1
MDL: MFCD00011076
Formula: AsF₆Li
Formula Weight: 195.85 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Gravimetric Analysis % AsF ₆ With Tetraphenylarsonium Chloride	94.0 - 98.9 %
X-Ray Diffraction	Conforms to Structure

Specification: PRD.0.ZQ5.10000029463

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.