

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
(3-Aminopropyl)triethoxysilane - 99%

Product Number: 440140
CAS Number: 919-30-2
MDL: MFCD00008207
Formula: C₉H₂₃NO₃Si
Formula Weight: 221.37 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared Spectrum	Conforms to Structure
Color Test	≤ 25 APHA
Purity (GC)	≥ 98.5 %
Expiration Date Period	-----
2 Years	

Specification: PRD.1.ZQ5.10000027944

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.