

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum - gauze, 100 mesh, 99.9% trace metals basis

Product Number: 298093 CAS Number: 7440-06-4 MFCD00011179 MDL:

Formula: Pt

Formula Weight: 195.08 g/mol

TEST Specification

Appearance (Color) Conforms

Silver

Appearance (Form) Conforms

Wire Mesh

0.076 mm

Diameter Conforms

ICP Major Analysis

Confirmed

Confirms Platinum Component

Purity Meets Requirements

99.9% Based On Trace Metals Analysis

Trace Metal Analysis < 1500.0 ppm

Specification: PRD.1.ZQ5.10000015472

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.