

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Tin(II) acetate

Product Number:

345164

CAS Number:

638-39-1

MDL:


MFCD00054408

Formula:

C₄H₆O₄Sn

Formula Weight:

236.8 g/mol



TEST

Specification

Appearance (Color)

White to Off-White, Tan or Yellow-Tan

Conforms to Requirements

Appearance (Form)

Powder and/or Chunks

Conforms to Requirements

Infrared spectrum

Conforms to Structure

Gravimetric Analysis

% Sn by Gravimetric Assay After Ignition to Oxide

49.0- 53.0 %

Remarks:

Specification Date : 11/29/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.