

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Boron triiodide - 95%

Product Number: 307629
CAS Number: 13517-10-7
MDL: MFCD00036286
Formula: BI₃
Formula Weight: 391.52 g/mol

BI₃

TEST	Specification
Appearance (Color)	White to Pink
Appearance (Form)	Powder or Crystals
Titration by AgNO ₃ % I	91.9 - 102.6 %

Specification: PRD.0.ZQ5.10000032334

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.