


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Methyl-3-nitrophenylboronic acid

Product Number: 521477
CAS Number: 80500-27-2
MDL: MFCD00191550
Formula: C₇H₈BNO₄
Formula Weight: 180.95 g/mol



TEST

Appearance (Color)
Appearance (Form)
Proton NMR spectrum
Purity (Titration by NaOH)

Specifi

White to
Powder or solid or crystals
Conforms to Structure
≥ 95.0 %

Specification: PRD.1.ZQ5.10000012937

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.