


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Gluconolactone - meets USP testing specifications

Product Number: G2164
CAS Number: 90-80-2
MDL: MFCD00006647
Formula: C₆H₁₀O₆
Formula Weight: 178.14 g/mol



TEST	Specification
Identity	Pass
Lead (Pb) < or = 0.001%	Conforms
Reducing Substances	Pass
Residual Solvent	Meets Requirements
Assay	99.0 - 101.0 %

Specification: PRD.1.ZQ5.10000044130

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.