

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification


Product Name:
N-(3,5-Bis(trifluoromethyl)phenyl)-2-chloroacetamide - ≥95%

Product Number: 912654

CAS Number: 790-75-0

Formula: C₁₀H₆ClF₆NO

Formula Weight: 305.60 g/mol



TEST	Specification
Appearance (Color)	Brown
Appearance (Form)	Powder
¹ H NMR Spectrum	Conforms to Structure
¹⁹ F NMR Spectrum	Conforms to Structure
Titration with AgNO ₃ After Oxygen Combustion	95.0 - 105.0 %
Purity (HPLC)	≥ 95.0 %

Specification: PRD.0.ZQ5.10000102733

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

