

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-Methylumbelliferyl phosphate (4-MUP) Liquid Substrate System - liquid

Product Number: **M3168**
CAS Number: 3368-04-5
MDL: MFCD00016969

Formula Weight: 256.15 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Solution
pH	9.8 - 10.2
Concentration	≥ 0.5 mM

Specification: PRD.0.ZQ5.10000050577

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.