


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N-Acetyl-L-cysteine - cell culture tested, BioReagent

Product Number: A9165
CAS Number: 616-91-1
MDL: MFCD00004880
Formula: C₅H₉NO₃S
Formula Weight: 163.19 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H ₂ O	
Carbon	36.0 - 38.0 %
Nitrogen	8.3 - 8.9 %
Infrared spectrum	Conforms to Structure
Purity (TLC)	≥ 99 %
Cell Culture Test	Pass

Specification: PRD.0.ZQ5.10000046255

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.