

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Retinyl acetate – synthetic, matrix dispersion, 475,000–650,000 USP units/g

Product Number:

R3250

CAS Number:

127-47-9

MDL:

MFCD00019413

Formula:


C₂₂H₃₂O₂

Formula Weight:

328.49 g/mol

Storage Temperature:

-20 °C



TEST	Specification
Appearance (Color)	Yellow
Appearance (Form)	Powder
Activity	475000 - 650000
USP units/gram based on HPLC comparison USP standard	

Specification: PRD.0.ZQ5.10000015190

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.