


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Pyruvic acid - 95%

Product Number: P76209
CAS Number: 127-17-3
MDL: MFCD00002585
Formula: C₃H₄O₃
Formula Weight: 88.06 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Purity (Titration by NaOH)
Recommended Retest Period
2 Years

Colorless to Light Yellow
Liquid
Conforms to Structure
94.0 - 106.0 %

Specification: PRD.0.ZQ5.10000018887

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.