

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium borohydride solution - SWS, ~12 wt. % in 14 M NaOH

Product Number:

452904



MDL:

MFCD00003518

Formula:

H4BNa

Formula Weight:

37.83 g/mol

TEST

Specification

Appearance (Color)

Colorless to Very Faint Yellow

Appearance (Form)

Viscous Liquid

Hydrogen Evolution

11.85 - 12.15 %

NaBH₄

Product of Supplier

Confirmed

Ascensus

Specification: PRD.2.ZQ5.10000011147

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.