


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Hydroxy-2-quinoxalinecarboxylic acid - 97%

Product Number: 108340
CAS Number: 1204-75-7
MDL: MFCD00006726
Formula: C₉H₆N₂O₃
Formula Weight: 190.16 g/mol



TEST	Specification
Appearance (Color)	Yellow to Gold
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Proton NMR Spectrum	Conforms to Structure
Purity (Titration by NaOH)	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000040903

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.