

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anti-UCP-3 antibody produced in rabbit – affinity isolated antibody, buffered aqueous solution

Product Number: U7757

TEST	Specification
Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	1.0 - 1.5 mg/ml
WB-Cell Line/Tissue Extract	293T Cells Expressing Human UCP3
WB-Western Blot Titer	1:1000-1:2000
IHC -Tissue	Human Tongue
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC-Titer	1:3000-1:6000
IHC -Tissue	Mouse Leg Muscle
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC-Titer	1:3000-1:6000
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000018087

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.