

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium hexafluoroaluminate - 99.98% trace metals basis

Product Number: 230049
CAS Number: 13775-53-6
MDL: MFCD00003507
Formula: AlF₆Na₃
Formula Weight: 209.94 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Chunks
Complexiometric EDTA (% Al)	12.0 - 14.0 %
X-Ray Diffraction	Conforms to Structure
ICP: Confirms Na and Al Compon	Confirmed
Purity	Meets Requirements
99.98% Based on Trace Metals Analysis	
Trace Metal Analysis	≤ 250.0 ppm

Specification: PRD.0.ZQ5.10000042240

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.